

Therapeutic drug monitoring (TDM) of thiopurine drugs

Kamisha L. Johnson-Davis, Ph.D., DABCC (CC,TC)

Associate Professor
Department of Pathology, University of Utah

Medical Director – Clinical Toxicology
ARUP laboratories




Topics

- Discuss the function of thiopurine drugs
- Thiopurine Drug Metabolism
- Pharmacogenetics of Thiopurine methyltransferase (TPMT)
- Describe analytical methods to support therapeutic drug management


Thiopurine drugs

Immunosuppressive antimetabolites


**Azathioprine
(Imuran®)**



**6-Mercaptopurine
(Purinethol®)**



**Thioguanine
(Tabloid®)**



↑
Analog of 6-Mercaptopurine

↙ ↘
These compounds are analogs of hypoxanthine and guanine

Thiopurine Drugs - Clinical Use

Acute lymphoblastic leukemia –

- Cancer in WBC
 - Most common cancer in children
 - 3,000 new cases / year (3.7-4.9 cases / 100,000) in children
 - 98% remission after treatment
-
- Reference: www.stjude.org

Thiopurine Drugs - Clinical Use

Inflammatory Bowel Disease


- Autoimmune Disease
 - 1.6 million Americans (80,000 children)
 - 70,000 new cases / year
 - Maintain remission 70-75%
-
- Reference: www.ccfa.org

Thiopurine Drugs - Clinical Use

Rheumatoid Arthritis

- Autoimmune Disease – tissues near joints
 - 1.3 million adults (>18 yr)
 - Incidence – 89 / 100,000 (65 – 74 yr);
41 / 100,000 (30 – 65 yr)
 - 50 – 60% remission after treatment
-
- Reference: www.cdc.gov




Metabolism of Thiopurine Drugs




Kennedy et al. 2013 Intern Med J. 43(3):278-286.

Thiopurine Drugs - cytotoxicity


- 6-TGN resembles DNA nucleotides
- Incorporates into replicating DNA
 - Blocks proliferation (T-cells, B-cells)
- Inhibits purine de novo synthesis (6-MMPN)
- Promote immunosuppression
- Induce T-cell apoptosis



6-TGMP



6-TGDP



6-TGTP

Is Thiopurine Therapy Effective in All Patients?

Table 2. Azathioprine tolerance: estimated proportions of patients

Please estimate the percentage of your patient population in which azathioprine is tolerated well, reasonably well and poorly.

Reported tolerance	Derma (n = 95) mean% (SD)	Gastro (n = 92) mean% (SD)	Rheum (n = 96) mean% (SD)	All (n = 283) mean% (SD)	P* (significant)
Well	69 (21)	75 (12)	59 (19)	67 (19)	<0.001 (yes)
Reasonably well	17 (17)	13 (9)	22 (12)	17 (14)	<0.001 (yes)
Poorly	14 (10)	13 (7)	19 (11)	15 (10)	<0.001 (yes)
Missing responses	4	1	6	11	

*ANOVA; P < 0.05.

Table 3. Estimated proportions of patients experiencing azathioprine side-effects

Which of the following side-effects have you had experience of in your patient population? Please estimate the percentage.



Side-effect	Derma (n = 95) mean% (SD)	Gastro (n = 92) mean% (SD)	Rheum (n = 96) mean% (SD)	All (n = 283) mean% (SD)	P* (significant)
Nausea	16 (16)	11 (10)	20 (18)	16 (15)	<0.001 (yes)
Pancreatitis	0.1 (0.5)	2 (2)	0.3 (1.3)	0.6 (1.7)	<0.001 (yes)
Neutropaenia	4 (6)	5 (4)	8 (11)	6 (8)	0.003 (yes)
Abnormal LFTs	9 (8)	5 (4)	10 (11)	8 (9)	<0.001 (yes)
Missing responses	2	8	7	17	

*ANOVA; P < 0.05.


Fargher et al. (2007) J. Clin Pharm Ther 32: 187-195

Thiopurine S-methyltransferase (TPMT)

- Cytoplasmic s-methylation of thiopurines
 - 28 kDa, soluble protein; 245 amino acids
- Sites of Metabolism
 - Heart, Liver, Kidney, Intestines, Pancreas
 - RBC and WBC
- Exhibits genetic polymorphisms
 - May lead to adverse effects from thiopurine drug use




Distribution of TPMT Activity – Caucasian Population




Shaeffeler E et al, Pharmacogenetics 2004


Metabolism of Thiopurine Drugs



Metabolism of Thiopurine Drugs



Metabolism of Thiopurine Drugs



TPMT Allele

Normal TPMT Allele



- TPMT is ~30 kb long
- Consists of 10 exons
 - 8 of which encode for protein
- Patients with 2 normal TPMT alleles
 - normal high activity


Shaeffeler E et al. 2004, Pharmacogenetics

Genetic Mutations in TPMT Alleles

- Single Nucleotide Polymorphisms (SNPs)
 - Occur in the coding region of the gene
 - Occur at exon-intron splice sites
 - Can decrease the function of the enzyme

- Variable number tandem repeats (VNTRs)
 - Occur in the promotor region of the gene
 - Decrease function of the enzyme
 - Most common – VNTR*3, *4, *5
 - Spire-Vayron et al, 1999; Roberts et al., 2008

Most Common Nonfunctional TPMT Alleles




- The four most common nonfunctional TPMT alleles account for >95% of low TPMT activity
- The most common allele variant found in the Caucasian population is the *TPMT* *3A allele
- The most common allele variant found in the Asian and African population is the *TPMT* *3C allele

Shaeffeler E et al. 2004, Pharmacogenetics

Mutations affect TPMT Activity

- Protein becomes unstable and undergo degradation
 - *TPMT*2*, **3A*, **3B*, **3C* have enhanced proteolysis of variant proteins
 - Tai et al., (1999) Pharmacogenetics 9:641-50.
- Decreased enzyme half-life
 - *TPMT*1* – **18 h** vs. *TPMT*2* and *TPMT*3* – **15 min**
 - Tai et al., (1997) Natl Acad Sci USA 94:6444-6449.
- Variant alleles cause decreased enzyme activity
 - Krynetski et al., (1995) Proc. Natl. Acad. Sci. 92:949-953.

Dose Adjustments for TPMT Activity



Relling MV et al. 2011 Clin Pharmacol Ther. 89(3):387-391.

Analytical Methods to support testing for thiopurine drugs

Testing Strategies to Access TPMT Activity

Phenotype

- Enzyme reaction – assess inactivation of thiopurine drug
- Quantification of thiopurine drug metabolites

Genotype


- Determines the allele pattern of a patient
- Primarily tests for 4 variant TPMT alleles –
 - *TPMT*2*, **3A*, **3B*, **3C*

Phenotype Assay – TPMT Enzyme


6-MP + SAM

TPMT

6-MMP



centrifuge




Plasma

White blood cells
& Platelets

Red blood cells

- Remove plasma & WBC
- Lyse RBCs with water



TPMT


- **Add cosubstrates:**
6-mercaptopurine & S-adenosylmethionine
- **Incubate reaction @ 37 °C for 1 h**
- **Stop reaction with acid or heat:**
- **Prepare sample for HPLC or LC/MS/MS**



HPLC



Triple-Quad MS



U/mL = nanomoles 6MMP

Result Interpretation

Low (< 17 U/mL) – high risk of bone marrow toxicity with standard dosing regimens; dose reduction of 80-90% is recommended

Intermediate (17.0 - 23.9) – risk of bone marrow toxicity with standard dosing regimens; dose reduction by 20-50% may be necessary

Normal (24 - 44 U/mL) – low risk of bone marrow toxicity with standard dosing regimens

High (>44 U/mL) – risk of therapeutic failure due to excessive inactivation of these drugs. Patients may require a higher dose


Phenotype test - Thiopurine Drug Metabolites

- Test quantifies 6-TGN and 6-MMPN levels in whole blood

Normal Reference Ranges

6-TGN: 230-400 pmol/ 8×10^8 red blood cells

6-MMPN: <5700 pmol/ 8×10^8 red blood cells



Phenotype Test

Pros

- Less expensive than genotype test
- May be more informative than Genotype

Cons

- Enzyme function can be altered by:
 - NSAIDS can inhibit enzyme activity
 - Specimen mishandling can impact TPMT enzyme activity
 - Blood transfusion
 - Not accurately assessed in patients who receive RBC transfusion within 60 days of testing


Genotype TPMT assay

- PCR
 - Mini-sequencing assay
 - Fluorescence detection
- Detect the most common alleles
 - *TPMT**2 [238G>C]
 - *TPMT**3A [460G>A and 719A>G]
 - *TPMT**3B [460G>A]
 - *TPMT**3C [719A>G]

Genetic Test - PCR

General OpenArray® Workflow

The samples and master mix are mixed on an Intermediate Sample Plate (up to half of plate)



Samples and 2x Master Mix loaded onto the OpenArray® plate.



Lid affixed onto array case, filled with fluid, sealed, and loaded onto carrier.

Plates are cycled and imaged on QS. (Optional: 9700 flat block for genotyping off-line cycling.)



Analyze



TaqMan Genotyper Software

ThermoFisher SCIENTIFIC

Genotype Test

Pros

- Genotype test is not affected by –
 - Blood transfusions
 - Concomitant medications

Cons

- Assay cannot detect all of the nonfunctional TPMT alleles (>30)
- Wildtype or Normal alleles are inferred if none of the four non-functional alleles are detected
- Cannot distinguish between a heterozygote or a compound heterozygote
 - *1/*3A or *3A/*3C, *3A/*3B, *3B/*3C

Summary

- Phenotype assay may provide more clinical information
 - Detects rapid metabolizer phenotypes and impaired enzyme activity from rare genetic variants
 - Genotype test may not detect all rare mutations
- Genetic testing can be performed before and after thiopurine drug therapy
 - Consider testing if erythrocyte TPMT activity is abnormal

Take Home Message

- Test patients for TPMT status prior to therapy
 - Phenotype
 - Genotype
- Patients with low or absent TPMT activity
 - Avoid Thiopurine therapy or reduce the standard dose
- Genotype/phenotype tests to assess TPMT activity do not replace routine clinical monitoring



Department of Pathology

© 2016 ARUP Laboratories