


# Applications for AI in the Clinical Parasitology Lab



BLAINE A. MATHISON, BS, M(ASCP)

TECHNICAL DIRECTOR OF PARASITOLOGY, ARUP LABORATORIES, SALT LAKE CITY, UT

ADJUNCT INSTRUCTOR, DEPT. OF PATHOLOGY, UNIVERSITY OF UTAH, SALT LAKE CITY, UT

# Disclosures in the Past Year

- Disclosures Related to this Talk:
  - » Techcyte Inc. (collaborator)
- Other Professional Disclosures:
  - » American Society for Microbiology (Editorial Board, *Journal of Clinical Microbiology*; Editorial Board, *ASM Case Reports*; honoraria and travel support)
  - » College of American Pathologists (Microbiology Committee Member, 2024-present)
  - » Illinois Society for Microbiology (travel support; honorarium)
  - » University of Minnesota (travel support)
  - » Association of Public Health Labs (travel support)
  - » Mayo Clinic (travel support)
  - » Southwestern Association for Clinical Microbiology (travel support)
  - » Elsevier (publication honorarium)

# Product Disclaimer

The mentioning or showing of commercial products in this talk does not constitute an endorsement nor renouncement of the products

# Learning Objectives

---

Review the Morphologic Identification of Intestinal Parasites

---

Describe Potential Workflow Changes when Adopting AI for Parasitology

---

Discuss Proficiency, Competency, and Quality Control as it Pertains to AI

---

# Review: Laboratory Diagnosis of Intestinal Parasites


---









# A Brief History of Laboratory Identification of Intestinal Parasites

- Microscopy
- Antigen Detection
- Molecular Detection
- Artificial Intelligence







# Review: Morphologic Identification of Intestinal Parasites

- Wet Mount (direct or concentrate)
- Trichrome/Iron Hematoxylin Smear
- Modified Acid-fast/Safranin
- Histopathology
- Gross Examination



# Review: Morphologic Identification of Intestinal Parasites

- **Wet Mount (direct or concentrate)**
- **Trichrome/Iron Hematoxylin Smear**
- **Modified Acid-fast/Safranin**
- Histopathology
- Gross Examination





# AI: Benefits, Challenges, and Practical Application

---



# Basic Concepts of the Application of AI

- First and most important, this is a screening tool, not a definitive parasite identification tool.
  - » Models were designed for sensitivity over specificity
- The application of the AI does not diminish the need for parasitology expertise, competency, and proficiency
- An experienced technologist must still evaluate images
- The goal is to reduce the number of manual reads for likely negative specimens
  - » **'negative screening tool'**

# Potential Benefits to Adopting AI


- Less time spent manually reading negative specimens
  - » In our lab, the negative percent for O&P is 92-95%
- More time devoted to scrutinizing positive or putative positive specimens
  - » Improves diagnostic skills and competence and confidence in parasite identification
- Increased sensitivity
- Less physical stressors and ergonomic issues
  - » Eye strain (microscopy), hand strain (microscope adjustment), poor posture
- Possibly more engaging for an increasingly younger workforce
- Potential for remote analysis
  - » CAUTION: could be constraints on where clinical specimens can be formally evaluated

# Workflow Challenges to Consider

- Slide preparation for both trichrome/MAF and wet mount
  - » Hills and valleys vs. homogenous smears (permanent smears)
  - » Coverslipping of permanent smears
  - » Wet mount preps to prevent drying and movement
- Labeling slides
- Choosing a scanner for volume and workflow
- Acquisition of specimens for validation
- Possible reorganization of workforce
  - » Technicians vs. technologists
- Quality control (upcoming section)
- Having back-up plans when there is failure at any step
  - » Scanner malfunctions
  - » Coverslipper breaks
  - » AI/upload issues


# Changes to Permanent Smear Processes

- Permanent smears (trichrome, MAF, SAF) will have to be prepared with a thin monolayer
  - » No more 'hills and valleys'
- Slides will need to be coverslipped prior to scanning
  - » Automated vs. manual coverslipping
- Slides will need a label with a barcode/QR code for reading by instrument




# Changes to Wet Mount Preps

- Precise measurements of stool and mounting medium (iodine solution) are required
  - » 20  $\mu$ L 50% Lugol's iodine + 50% 10% glycerol in PBS
  - » 15  $\mu$ L stool specimen
- Use of a template to make sure the stool is in the designated scan area
- Slides will need a bar code/QR code for reading
- The preparation of smears will have to be carefully timed to prevent drying
- Make sure there are no air bubbles or parts of the coverslip devoid of stool/mounting medium
- Challenging specimens may not be scannable and will have to be manually read per normal
  - » Grainy, excess mucus, 'interfering substances'





# Selecting a Scanner



Pramana SpectralHT2



Pramana SpectralM-Pro



Hamamatsu  
NanoZoomer 360



Grundium Ocus40

**NOTE:** This slide is not an endorsement or renouncement of shown products

# Quality Control, Proficiency, Competency

---




# Competency and Proficiency

- Competency = tests the individual's skill set
  - » In our lab, an employee must be competent on traditional microscopic identification of the parasites prior to learning AI
  - » After initial training new employees get unknowns consisting of a mix of positive and negative scans
    - If organisms are missed, scans are reviewed with Education Coordinator
    - If several organisms are missed, a second (different) test is given
  - » In our lab, we have quarterly competency testing both traditional parasite ID and recognition of parasites via AI
- Proficiency = tests the workflow
  - » Technically, proficiency specimens should be handled in the same manner of clinical specimens
  - » For proficiency specimens to be scanned, scan area must be correlated with premade proficiency specimens (e.g., trichrome)

# Quality Control – Permanent Smears

- QC is performed daily on trichome slides.
  - » QC slide for scanning is not the same for the stain itself
- For the scanning QC slides, both positive and negative slides are processed; a lot for the positive smear is made from a specimen containing at least one protozoan.
- Ranges are established with each **slide** (not lot) and counted boxes (identified by the software)
  - » If total counted boxes for a given organism are outside of  $\pm 3$  SD range, the slide is rescanned
  - » If after rescanning the total counted boxes are outside of the  $\pm 3$  SD range, the QC is acceptable **if** the first 50 boxed organisms are correct **and** >80% of the prevalence regions are valid (in focus)
- Slides are retired after six months to prevent bleaching from scanner light

# Wet Mounts





# Case Examples – How does it Work?

---





# What can the AI detect? Permanent Smears

## Trichrome

- » *Giardia duodenalis* – cysts and trophozoites
- » *Entamoeba* sp. non-*hartmanni* - trophozoites
- » *Entamoeba hartmanni* - trophozoites
- » *Dientamoeba fragilis*, *Iodamoeba buetschlii*, *E. nana* – trophozoites
- » *Blastocystis* species
- » *Chilomastix mesnili* – trophozoites
- » *Cyclospora* species – trophozoites
- » White Blood Cells
- » Red Blood Cells

## Modified Acid-Fast

- » *Cryptosporidium* species
- » *Cyclospora* species

Yeast was trained as an anti-class and available for those that report yeast in OPs

# What can the AI detect? Wet Mounts

## Protozoans


- » *Giardia duodenalis* – cysts and trophozoites
- » *Entamoeba* sp. non-*hartmanni* – cysts and trophozoites
- » *Blastocystis* species
- » *Chilomastix mesnili* – cysts
- » *Cyclospora* species – oocysts
- » *Cystoisospora belli* - oocysts
- » *Balantioides coli* – cysts and trophozoites
- » *Endolimax nana* – cysts
- » *Iodamoeba buetschlii* – cysts
- » Misc. Small Protozoans\*

\**Iodamoeba buetschlii*, *Dientamoeba fragilis*, *Endolimax nana*, *Chilomastix mesnili*, and *Entamoeba hartmanni* trophozoites; *E. hartmanni* cysts

## Helminths

- » *Enterobius vermicularis*
- » *Ascaris lumbricoides*
- » *Trichuris trichiura*
- » Hookworm/*Trichostrongylus*
- » *Strongyloides* species – L1 larvae
- » *Paracapillaria philippinensis*
- » *Taenia* species
- » Fish Tapeworm (*Diphyllbothrium*-complex)
- » *Rodentolepis nana*
- » *Hymenolepis diminuta*
- » *Schistosoma mansoni*
- » *Schistosoma japonicum*
- » *Fasciola hepatica*/*Fasciolopsis buski*
- » *Clonorchis*/*Opisthorchis* species
- » *Paragonimus* species

# Using the AI interface



Search

BL Blaine Mathison (Technologist Role)  
ARUP Laboratories - PAFT

RUO

Worklist +

2.1. OPFEC Pending Needs Assignment

139 512

1. Assigned to Me 2.1. OPFEC Pending Needs Assignment 2.2. PARAST Pending Needs Assignment 3.1. QC ready 6. Duplicate scans 7. Reviewed by Me

☐ paft\_op\_24472

☐

Human Fecal Ova Parasite 10/29/2025

0 of 2 scans reviewed

Pending review  
Unassigned

☐

Human Fecal Ova Parasite 10/29/2025

0 of 2 scans reviewed

Pending review  
Unassigned

☐

Human Fecal Ova Parasite 10/29/2025

0 of 2 scans reviewed

Pending review  
Unassigned

☐

Human Fecal Ova Parasite 10/29/2025

0 of 2 scans reviewed

Pending review  
Unassigned

☐

Human Fecal Ova Parasite 10/29/2025

0 of 2 scans reviewed

Pending review  
Unassigned

# Using the Interface





25087149329-OPFEC

Human Fecal Ova Parasite





## Case details

Case ID	Case number	Case status	Assigned to
[REDACTED]	[REDACTED]	Reviewed	Technologist [REDACTED]
Created at	Completed at	Archived at	Origin
04/02/2025, 6:46:50 AM	04/07/2025, 6:49:41 AM	Not set	External
Result	Patient name	Sex	Date of birth
Positive			
Batch ID	Shelf	Quality control	
paft_op_21709	Archived	False	


## Trichrome

Scan label	Sample type	Barcode	Status	Actions
[REDACTED]	Human Fecal	[REDACTED]	Reviewed	   

## Wet Mount

Scan label	Sample type	Barcode	Status	Actions
[REDACTED]	Human Fecal W...	[REDACTED]	Reviewed	   

# Trichrome smear



Search

BL Blaine Mathison (Technologist Role)  
ARUP Laboratories - PAFT

Barcode [redacted] A\_scanner Case type Human Fecal Ova Parasite Sample type Human Fecal

Fecal Smear

381

Blastocystis sp.

Prevalence \*

1+ 142 ☒

Giardia duodenalis

8 ☐

D. fragilis/E. nana/I. buetschlii

80 ☐

Entamoeba hartmanni

48 ☐

Entamoeba sp.

35 ☒













Chilomastix mesnili

38 ☐


Cyclospora sp.


12 ☐

Blastocystis sp.



Scan list

Trichrome 1 total 

Wet Mount 1 total 

Assessment


Cannot edit scan Assessment already finished

Previous assessment Technologist [redacted] days ago)


Notes

ARUP LABORATORIES | UNIVERSITY OF UTAH HEALTH | 25

Example: *Blastocystis* sp.




## Examples





# Trichrome smear



Search

BL Blaine Mathison (Technologist Role)  
ARUP Laboratories - PAFT

Barcode [REDACTED] A\_scanner Case type Human Fecal Ova Parasite Sample type Human Fecal

D. fragilis/E. nana/I. buetschlii  
80 ☐

Entamoeba hartmanni  
48 ☐

Entamoeba sp.  
35 ☒

Chilomastix mesnili  
38 ☐


Cyclospora sp.  
12 ☐

White blood cells  
2 ☐


Red blood cells  
16 ☐


Prevalence Region OPFEC  
11

Entamoeba sp.



Scan list

Trichrome  
1 total  


Wet Mount  
1 total  


Assessment


Cannot edit scan Assessment already finished

Previous assessment Technologist [REDACTED] (206 days ago)

Notes

ARUP LABORATORIES | UNIVERSITY OF UTAH HEALTH | 27

# Example: *Entamoeba* sp.





# Trichrome – Prevalence Regions

**D. fragilis/E. nana/I. buetschlii** 80 ☐

**Entamoeba hartmanni** 48 ☐

**Entamoeba sp.** 35 ☒

**Chilomastix mesnili** 38 ☐


**Cyclospora sp.** 12 ☐

**White blood cells** 2 ☐


**Red blood cells** 16 ☐


**Prevalence Region OPFEC** 11

Prevalence Region OPFEC



Scan list
Assessment

Trichrome 1 total 

Wet Mount 1 total 

**Cannot edit scan** Assessment already finished

Previous assessment *Technologist* (206 days ago)

Notes

# Slide View

D. fragilis/E. nana/I. buetschlii80☐

Entamoeba hartmanni48☐

Entamoeba sp.35☒

Chilomastix mesnili38☐


Cyclospora sp.12☐

White blood cells2☐

Red blood cells16☐

Prevalence Region OPFEC11

Entamoeba sp.



Scan list

Trichrome1 total☒

Wet Mount1 total☒


Assessment

Cannot edit scan Assessment already finished

Previous assessment *Technologis* (206 days ago)

Notes

10x 20x 40x 50x 80x 100x



2 mm 0.58x


Read only

ARUP LABORATORIES | UNIVERSITY OF UTAH HEALTH | 30



# Slide View

Entamoeba sp.




20 µm 82x

Scan list Assessment

Microscopic image showing Entamoeba sp. (a large, pear-shaped trophozoite with a prominent nucleus and a smaller, darker nucleus) and other cellular structures. The image is displayed in a slide view interface with a toolbar at the top and a navigation bar at the bottom.

# Wet Mount



Search

BL Blaine Mathison (Technologist Role)  
ARUP Laboratories - PAFT

Barcode [REDACTED] Slot A-06 Cassette UNK Case type Human Fecal Ova Parasite Sample type Human Fecal Wet Mount

Human Fecal Wet Mount

842


Blastocystis spp.  
10 ☒

Cyclospora spp.  
4 ☐

Ascaris lumbricoides  
828 ☒

Prevalence Region HFW  
14

Blastocystis spp.



Scan list

Trichrome [REDACTED]  
1 total ☒

Wet Mount [REDACTED]  
1 total ☒

Assessment

Cannot edit scan Assessment already finished


Previous assessment Technologist [REDACTED] 206 days ago

Notes


Read only



# Example: *Blastocystis* sp.



# Wet mount



Search

BL Blaine Mathison (Technologist Role)  
ARUP Laboratories - PAFT

S1 Slot A-06 Cassette UNK Case type Human Fecal Ova Parasite Sample type Human Fecal Wet Mount

Human Fecal Wet Mount

842


Blastocystis spp.  
10 ☒

Cyclospora spp.  
4 ☐


Ascaris lumbricoides  
828 ☒


Prevalence Region HFW  
14

Ascaris lumbricoides



Scan list

Trichrome  
1 total 

Wet Mount  
1 total 

Assessment


Cannot edit scan Assessment already finished

Previous assessment *Technologist* 206 days ago

Notes

Read only


# Example: *Ascaris lumbricoides*






# Some fun examples since going live with wet mount:


*Schistosoma mansoni*




*Taenia* spp.



*Strongyloides* spp.




*Trichuris trichiura*



*Enterobius vermicularis*



*Rodentolepis nana*



# A clean negative...

Case type Human Fecal Ova Parasite Sample type Human Fecal

Fecal Smear 149

Blastocystis sp. 18 ☐

Giardia duodenalis 22 ☐

D. fragilis/E. nana/I. buetschlii 3 ☐

Entamoeba sp. 2 ☐

Chilomastix mesnili 6 ☐

Cyclospora sp. 14 ☐

White blood cells 1 ☐

Blastocystis sp.

Giardia duodenalis

Trichrome 1 total

Wet Mount 1 total

Assessment (1 of 521)

Cannot edit scan Not assigned to you

Notes

Save



# And its corresponding wet mount...

techcyte

Search


BL Blaine Mathison (Technologist Role) ARUP Laboratories - PAFT

Slot A-01 Cassette UNK Case type Human Fecal Ova Parasite Sample type Human Fecal Wet Mount


human Fecal Wet Mount

Blastocystis spp.	17	
Giardia duodenalis	2	
Cyclospora spp.	1	
Entamoeba species	1	
Misc. Small Protozoans	1	
Prevalence Region HFW	14	


Blastocystis spp.



Giardia duodenalis



Cyclospora spp.



Scan list

Trichrome	1 total
Wet Mount	1 total

Assessment ( 1 of 521 )


Cannot edit scan Not assigned to you

Notes

Save



# Kinyoun's Modified Acid-fast Stain



BL Blaine Mathison (Technologist Role)  
ARUP Laboratories - PAFT

B\_scanner Case type Human Fecal Modified Acid Fast Sample type Human Fecal Modified Acid Fast (PARAST)

count total present

Human Fecal Modified Acid Fast (PARAST)

0 40


Cryptosporidium sp. oocysts

40 40 ☒

Prevalence Region OPFEC

11 11

### Cryptosporidium sp. oocysts



Assessment


Cannot edit scan Assessment already finished

Previous assessment Technologist (6 days ago)


Notes

Read only

# Example: *Cryptosporidium* sp.




# Lab Workflow



# Would the lab have missed this?

	count	total	present
Fecal Smear	0	42	
Blastocystis sp.	2		<input type="checkbox"/>
Giardia duodenalis	0	1	
Giardia Troph	1		<input type="checkbox"/>
Entamoeba sp.	0	2	
Entamoeba Troph	2		<input type="checkbox"/>
White blood cells	4		<input type="checkbox"/>
Red blood cells	33		<input type="checkbox"/>
Prevalence Region OPFEC	11		


Blastocystis sp.




Giardia Troph



Entamoeba Troph



Giardia Troph



# Or what about this one?

Human Fecal Wet Mount

14

Giardia duodenalis

11

☐

Entamoeba species

1

☐

Iodamoeba buetschlii

1

☐

Trichuris trichiura


1

☐

Prevalence Region HFW

14

Giardia duodenalis



Entamoeba species

Assessment

Cannot edit scan Assessment already finished

Previous assessment *Technologist* (08 days ago)

Notes


Read only




# But wait, there's more!

Human Fecal Wet Mount	14	
Giardia duodenalis	11	<input type="checkbox"/>
Entamoeba species	1	<input type="checkbox"/>
Iodamoeba buetschlii	1	<input type="checkbox"/>
Trichuris trichiura	1	<input type="checkbox"/>
Prevalence Region HFW	14	

Trichuris trichiura



Prevalence Region HFW



Assessment

**Cannot edit scan** Assessment already finished


Previous assessment *Technologist* (redacted) 08 days ago)

Notes


# And another challenging case...

Human Fecal Wet Mount	8	
Blastocystis spp.	3	<input type="checkbox"/>
Cyclospora spp.	1	<input type="checkbox"/>
Entamoeba species	3	<input type="checkbox"/>
Strongyloides spp.	1	<input type="checkbox"/>
Prevalence Region HFW	14	

Strongyloides spp.




Prevalence Region HFW



Assessment

Cannot edit scan Not assigned to you





# In Closing...

---



# In Closing. AI can improve the workflow in a parasitology lab by....

- Improving turnaround time (reducing the time needed to manually read negative specimens)
- Screening out negative specimens
- Reduce ergonomic injuries associated with prolonged microscopy
- Improve employee satisfactions (spending more time reading and confirming positive slides)
- Increased competence and confidence in parasite identification
- More engaging for an increasingly younger workforce

# Acknowledgements

## ARUP – Medical Directorship/R&D/R&I

Dr. Marc R. Couturier  
Dr. Adam Barker  
Dr. Mark Astill  
Dr. Jeff Stevensen  
Weston Hymas  
Jessica Kohan  
Katie Knight

## ARUP – Parasitology Laboratory

Cole Neider  
Ryan Jensen  
Madison Sant  
Alicia Backman

## Techcyte, Inc.

Ben Cahoon  
Rick Smith  
John Walker  
Brian Cahoon  
Jill Potts  
Ben Black  
Falon Markow  
Amy Wood

Blaine Mathison [blaine.mathison@aruplab.com](mailto:blaine.mathison@aruplab.com)





*ARUP is a nonprofit enterprise of the University of Utah and its Department of Pathology.*